

DORO PHONE SECURE 780X IUP WHITE



Réf. HDSECU780P



DORO PHONE SECURE 780X IUP WHITE



Réf. HDSECU780P

FEATURES			
Medical device	No		
Marquage	CE	Barcode	7322460079826
Guarantee (month):	24	Customs Code	85171200
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	16,5
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,336
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



11/07/2022 21:12 - This article is available on our website : www.gsh-med.fr with the following reference HDSECU780P