

SMARTPHONE DORO 8100 PLUS, GRAPHITE COLOR



Réf. HDSMAR01



SMARTPHONE DORO 8100 PLUS, GRAPHITE COLOR



Réf. HDSMAR01

FEATURES			
Medical device	No		
Marquage	Barcode	7322460082321	
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,3	Width packaging (cm):	
Length (cm) :	15,6	Length packaging (cm) :	
Height (cm) :	0,9	Height packaging (cm) :	
Weight (kg) :	0,167	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



17/05/2022 12:57 - This article is available on our website : www.gsh-med.fr with the following reference HDSMAR01