



DORO WATCH BLACK



Réf. HDWATC01N



DORO WATCH BLACK



Réf. HDWATC01N

FEATURES			
Medical device	No		
Marquage		Barcode	7322460083595
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/12/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference HDWATC01N