

DORO COMFORT 1015 WHITE, WITH ANSWERING MACHINE



Réf. HDCOMF06W



DORO COMFORT 1015 WHITE, WITH ANSWERING MACHINE



Réf. HDCOMF06W

FEATURES			
Medical device	No		
Marquage	Barcode	7322460060473	
Guarantee (month):	24	Customs Code	8517110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	20,0	Width packaging (cm):	7,0
Length (cm) :	16,0	Length packaging (cm) :	16,0
Height (cm) :	6,9	Height packaging (cm) :	20,0
Weight (kg) :	0,500	Weight packaging (kg) :	0,550
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



25/08/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference HDCOMF06W