

DUO DORO COMFORT 1015 WHITE, WITH ANSWERING MACHINE



Réf. HDCOMF07W



DUO DORO COMFORT 1015 WHITE, WITH ANSWERING MACHINE



Réf. **HDCOMF07W**

FEATURES			
Medical device	No		
Marquage	Barcode	7322460060503	
Guarantee (month):	24	Customs Code	8517110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,0	Width packaging (cm):	
Length (cm) :	18,0	Length packaging (cm) :	
Height (cm) :	16,0	Height packaging (cm) :	
Weight (kg) :	0,850	Weight packaging (kg) :	0,846
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



05/10/2022 12:36 - This article is available on our website : www.gsh-med.fr with the following reference HDCOMF07W