


MOBI CORAL RED WITH ADULT CUFF (M)



Réf. 513701



CE
0459



MOBI CORAL RED WITH ADULT CUFF (M)



Réf. 513701

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446027964
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,520
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



09/12/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference 513701

All characteristics are given for information only and may be changed