



MOBI POWDERY PINK WITH ADULT CUFF (M)



Réf. 513703



CE
0459



MOBI POWDERY PINK WITH ADULT CUFF (M)



Réf. 513703

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446027988
Guarantee (month):	36	Customs Code	90189100
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,520
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



09/12/2022 12:39 - This article is available on our website : www.gsh-med.fr with the following reference 513703

All characteristics are given for information only and may be changed