

LIAN CUFF/1 TUBE/NYLON (WITH CONNECTOR)/BLACK/TPU BLADDER/ LARGE (L)

Réf. 334106



LIAN CUFF/1 TUBE/NYLON (WITH CONNECTOR)/BLACK/TPU BLADDER/ LARGE (L)



Réf. **334106**

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | | Barcode | 3700446026974 |
| Guarantee (month): | 6 | Customs Code | 9018901000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm) : | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,215 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



26/04/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference 334106

All characteristics are given for information only and may be changed