

LIAN CUFF/1 TUBE/NYLON (WITH CONNECTOR)/GREY/ADULT (M)

Réf. 334111



LIAN CUFF/1 TUBE/NYLON (WITH CONNECTOR)/GREY/ADULT (M)



Réf. 334111

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446025052
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:17 - This article is available on our website : www.gsh-med.fr with the following reference 334111

All characteristics are given for information only and may be changed