

KIT 1 TUBE S-XS-XXS CARTOON CUFFS WITH CONNECTOR/CHILD

Réf. 334150



KIT 1 TUBE S-XS-XXS CARTOON CUFFS WITH CONNECTOR/CHILD



Réf. 334150

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446026622
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,335
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 334150

All characteristics are given for information only and may be changed