

## LIAN CUFF/1 TUBE/COTTON (WITH CONNECTOR)/NAVY BLUE/LARGE ADULT (L)

Réf. 336402



## LIAN CUFF/1 TUBE/COTTON (WITH CONNECTOR)/NAVY BLUE/LARGE ADULT (L)



Réf. **336402**

Cotton cuff with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446028527
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,330
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 336402

All characteristics are given for information only and may be changed