

# X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W

Réf. OTA1002



# X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W



Réf. OTA1002

These negatoscopes feature an X-ray mounting system for easy installation and removal, without affecting the X-ray. They can be placed on a table or mounted on a wall.

- The screen structure is made of 10 cm thick, high-strength aluminum alloy.
- The corners are protected by an ABS-coated steel plate.
- Methacrylic polymer screen
- Negotioscope is in high-strength aluminum alloy
- Interior and exterior coated with antistatic spray powder
- Voltage / Frequency: 220V / 50-60Hz
- Luminance: 2,000 cd/m<sup>2</sup>
- Average illumination > 6300 lux
- Lighting life: 50,000 h
- Color temperature: 10,000 k
- Exterior dimensions: 80 x 56 x 10 cm
- Screen dimensions: 72 x 43.5 cm

FEATURES					
Medical device	No				
Marquage	CE	Barcode	3661040195865		
Guarantee (month):	12	Customs Code	90105010		
Sterile	No				
INFORMATIONS		Main unit		Packaged unit	
Width (cm) :		Width packaging (cm):		68,2	
Length (cm) :		Length packaging (cm) :		87,0	
Height (cm) :		Height packaging (cm) :		15,1	
Weight (kg) :	8,500	Weight packaging (kg) :		12,000	
SPECIFIC STORAGE CONDITION					
Temperature min	(N/A)	Temperature max		(N/A)	
Humidity min	(N/A)	Humidity max		(N/A)	



11/07/2023 16:51 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OTA1002

All characteristics are given for information only and may be changed