

## AUDIOMETER SIBELMED SIBELSOUND DUO A OPTION AUTO HUGHSON WESTLAKE

Réf. DSI009133



# AUDIOMETER SIBELMED SIBELSOUND DUO A OPTION AUTO HUGHSON WESTLAKE

Réf. **DSI009133**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0318	Barcode	
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	32,0	Width packaging (cm):	32,0
Length (cm) :	44,0	Length packaging (cm) :	44,0
Height (cm) :	14,0	Height packaging (cm) :	14,0
Weight (kg) :	2,700	Weight packaging (kg) :	2,700
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



29/06/2022 12:23 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSI009133