

## ARRAL TROUSERS, WHITE, SIZE 5, HOLTEX



Réf. WSARR07\_95



## ARRAL TROUSERS, WHITE, SIZE 5, HOLTEX



Réf. WSARR07\_95

FEATURES			
Medical device	No		
Marquage	Barcode	3661040198026	
Guarantee (month):	12	Customs Code	62034311095
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	31,0	Width packaging (cm):	31,0
Length (cm) :	40,0	Length packaging (cm) :	40,0
Height (cm) :	1,5	Height packaging (cm) :	1,5
Weight (kg) :	0,236	Weight packaging (kg) :	0,236
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



29/06/2022 21:05 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference WSARR07\_95