

ARRAL TROUSERS, WHITE, SIZE 6, HOLTEX



Réf. WSARR07_96



ARRAL TROUSERS, WHITE, SIZE 6, HOLTEX



Réf. WSARR07_96

FEATURES			
Medical device	No		
Marquage	Barcode	3661040198033	
Guarantee (month):	12	Customs Code	62034311095
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	31,0	Width packaging (cm):	31,0
Length (cm) :	40,0	Length packaging (cm) :	40,0
Height (cm) :	1,5	Height packaging (cm) :	1,5
Weight (kg) :	0,236	Weight packaging (kg) :	0,236
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:06 - This article is available on our website : www.gsh-med.fr with the following reference WSARR07_96