

LIMB ELECTRODES CLAMP FOR ECG CARDIOMATE X 4

Réf. 260021



LIMB ELECTRODES CLAMP FOR ECG CARDIOMATE X 4



Réf. 260021

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0123		Barcode	3700446025519
Guarantee (month):	60		Customs Code	9018110000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,200	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



21/11/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference 260021

All characteristics are given for information only and may be changed