

LIAN NANO CLINIC/SILVER RING / MULTI CUFFS (L, M, S)

Réf. 513270



CE
0459



LIAN NANO CLINIC/SILVER RING / MULTI CUFFS (L, M, S)



Réf. 513270

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024741
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,873
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:28 - This article is available on our website : www.gsh-med.fr with the following reference 513270

All characteristics are given for information only and may be changed