

LIAN NANO/SILVER RING/ CARTOON DESIGN COTTON / CHILD (S)

Réf. 513229



CE
0459



LIAN NANO/SILVER RING/ CARTOON DESIGN COTTON / CHILD (S)



Réf. 513229

Lightweight and robust hand-held, combines style and performance. Lightweight, shock-resistant ABS casing. Easy-read large dial. Cartoon-pattern cotton cuffs with screw-on metal quick coupling. Ergonomic inflation bulb and brass valve with ABS cap. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446007195
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,5
Length (cm) :		Length packaging (cm) :	19,5
Height (cm) :		Height packaging (cm) :	7,0
Weight (kg) :	0,500	Weight packaging (kg) :	0,526
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



26/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 513229

All characteristics are given for information only and may be changed