

LIAN NANO CLINIC/SILVER RING / LARGE ADULT (L)

Réf. 513250



LIAN NANO CLINIC/SILVER RING / LARGE ADULT (L)



Réf. 513250

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe I
Marquage			Barcode 3700446024710
Guarantee (month):	12	Customs Code	90189010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) : 0,565	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:27 - This article is available on our website : www.gsh-med.fr with the following reference 513250

All characteristics are given for information only and may be changed