

ALUMINIUM BOX 180 X 90 X 30 MM - GOLD



Réf. L30211



ALUMINIUM BOX 180 X 90 X 30 MM - GOLD



Réf. **L30211**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446000752
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,151
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/10/2022 12:43 - This article is available on our website : www.gsh-med.fr with the following reference L30211

All characteristics are given for information only and may be changed