

SUCCION UNIT NEW ASKIR 230/12V, 2000ML FLOVAC DISPOSABLE LINER

Réf. BCASPI23

CA-MI

GARANTIE
3
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
II
a



SUCCION UNIT NEW ASKIR 230/12V, 2000ML FLOVAC DISPOSABLE LINER

CA-MI

Réf. **BCASPI23**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139098
Guarantee (month):	36	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,5	Width packaging (cm):	
Length (cm) :	38,5	Length packaging (cm) :	
Height (cm) :	30,0	Height packaging (cm) :	
Weight (kg) :	4,650	Weight packaging (kg) :	5,000
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



07/05/2021 16:28 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI23

All characteristics are given for information only and may be changed

2/2