

LUER ST STEEL TRACHEAL TUBE 2 PCS N°2 DIAM. 4 MM

Ref. IC43200



26/08/2024 22:28:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

LUER ST STEEL TRACHEAL TUBE 2 PCS N°2 DIAM. 4 MM



Ref. IC43200

Features

CE Marking	CE	Medical device	Class IIa
Bar code	3661040083810	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	5.4
Packed unit / Weight	0.025		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 22:28:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr