

## NEW EMIVAC PORTABLE FOOT- OPERATED SUCTION UNIT




Ref. BCASPI75



 29/08/2024 15:08:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## NEW EMIVAC PORTABLE FOOT-OPERATED SUCTION UNIT



Ref. BCASPI75

**Pedal surgical Suction Unit mobile that is suitable for emergency aspirations. Requiring no lubrication or battery change, it is therefore always ready to use.**

### Description

Its use is suitable for gentle suctions such as tracheal aspiration as well as young children, up to 300 mmHg.

Main applications: Car-Ambulance, Emergency Department, Mobile Hospital, Obstetrics- Supplied with: an autoclavable

collection jar of 400 ml in polycarbonate and with an anti-reflux valve / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient)

- Motor: manual on foot

- Maximum vacuum (adjustable): 0.40 Bar -40 kPa -300 mmHg

- Max free air flow: 22 l/min

- Weight: 1.15 kg

- Dimensions: 22 x 16 x 8 cm

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8033622540175	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	16
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	22
Packed unit / Weight	1.6	Bare Unit / Weight (kg)	1.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:08:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)