

COLLYER PELVIMETER (GRADUATED IN CM AND INCHES)



Ref. IP45901



27/08/2024 04:23:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

COLLYER PELVIMETER (GRADUATED IN CM AND INCHES)



Ref. IP45901

Features

CE Marking	CE	Medical device	Class Im (measuring function)
Bar code	3661040026046	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	15
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	33.7
Packed unit / Weight	0.253	Bare Unit / Weight (kg)	0.253
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 04:23:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr