

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 200 MM - 2 MM



Ref. IP46844



27/08/2024 04:37:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 200 MM - 2 MM



Ref. IP46844

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040026442	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	12.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	29	Bare Unit / Length (cm)	29.2
Packed unit / Weight	0.207	Bare Unit / Weight (kg)	0.207
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 04:37:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr