

SAUERBRÜCK BONE RONGEUR MULTIACTION STRAIGHT CURVED JAW 30 CM



Ref. IP49430



27/08/2024 04:57:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SAUERBRÜCK BONE RONGEUR MULTIACTION STRAIGHT CURVED JAW 30 CM



Ref. IP49430

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040027388	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	17.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	31
Packed unit / Weight	0.643	Bare Unit / Weight (kg)	0.643
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 04:57:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr