

CATHETER FOR JAR DIAM 54X400MM NEW HOSPIVAC



Ref. BCASPI130



29/08/2024 17:37:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CATHETER FOR JAR DIAM 54X400MM NEW HOSPIVAC



Ref. BCASPI130

Description

- Polycarbonate pipe 54 mm, 400 mm long.
- Fully autoclavable (121°C – 15 min)

Features

CE Marking	CE	Medical device	No
Bar code	0000000000789	HS Code	4202121900
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	43	Packed unit / Weight	0.256
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:37:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr