

AUTOCLAVABLE 4000 / 5000ML JAR HOLDER NEW HOSPIVAC



Ref. BCASPI132



29/08/2024 17:39:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

AUTOCLAVABLE 4000 / 5000ML JAR HOLDER NEW HOSPIVAC



Ref. BCASPI132

Features

CE Marking	NONE	Medical device	No
Bar code	0000000000802	HS Code	4202121900
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	6
Packed unit / Length (cm)	40	Packed unit / Weight	0.388
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		


29/08/2024 17:39:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr