

STAINLESS STEEL SCREEN 2 PANELS WITH STRAINED WHITE CURTAINS



Ref. OS0062



 22/11/2024 20:48:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

STAINLESS STEEL SCREEN 2 PANELS WITH STRAINED WHITE CURTAINS



Ref. OS0062

Description

- Panels (155 x 61 cm) with fixed castors without brakes
- Stainless steel tube frame
- Washable plastic curtains
- Fire resistance class: M3
- 4 double castors. Diameter: 60 mm

Features

CE Marking	NONE	Medical device	No
Bar code	3661040162621	HS Code	9403899000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	64	Packed unit / Height (cm)	9
Packed unit / Length (cm)	177	Packed unit / Weight	7
Bare Unit / Weight (kg)	6		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

22/11/2024 20:48:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr