

INSTRUMENT TROLLEY 2 SHELVES WITH DRAWER 60X40CM

Ref. OSA6042



29/08/2024 18:27:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

INSTRUMENT TROLLEY 2 SHELVES WITH DRAWER 60X40CM



Ref. OSA6042

High quality pedestal table, ideal for all types of medical practices or hospital services

Description

- Round tube structure in non-magnetic stainless steel for a better grip
- Stainless steel trays
- Removable galleries- Silent rubber wheels
- Permissible load: 80 kg
- Wheel diameter: 75 mm with brakes
- Drawer under upper tray

Drawer dimension: Length x Width x Height
36 x 41.7 x 9 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168302	HS Code	9402900000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	45	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	21	Bare Unit / Height (cm)	80
Packed unit / Length (cm)	85	Bare Unit / Length (cm)	60
Packed unit / Weight	12.4	Bare Unit / Weight (kg)	10.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 18:27:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr