

LIAN METAL HAND-HELD MANOMETER



Ref. 310045



30/08/2024 17:20:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

LIAN METAL HAND-HELD MANOMETER



Ref. 310045

Description

Dials 0 to 300 mmHg

Features

CE Marking	NONE	Medical device	Class Im (measuring function)
Bar code	3700446000141	HS Code	90189010
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.38
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:20:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr