

OXYMETER FINGER TYPE WITH CHILD SENSOR "DEPORTED" VELCRO PC-60E



Ref. DPOXY03



29/08/2024 16:16:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

OXYMETER FINGER TYPE WITH CHILD SENSOR "DEPORTED" VELCRO PC- 60E



Ref. DPOXY03

Description

The PC-60E suitable for use in the hospital, medical office and home.

- Automatic activation upon insertion of the finger into the oximeter
- Automatic shutdown after 6 seconds without activity
- OLED display resolution of SpO2, PR, bargraphand plethyspogramme
- Low battery indicator Visual
- Color Display Multi
- Audible and visual alarms
- Suitable for infants with Velcro sensor
- Dimensions: 66 x 36 x 33mm
- Weight: 60g (including batteries)
- Battery life: 30 hours continuous use
- Supplied with 2 AAA batteries, wrist strap 1 and 1 storage bag

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6949006031103	HS Code	9018909911
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	6.1
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	3.3
Packed unit / Length (cm)	22	Bare Unit / Length (cm)	3.1
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.059
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


29/08/2024 16:16:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr