

## XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III

# KERN®

Ref. EKMWS07\_CIII



CE  
0297

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
3 ANS  
**3**  
GUARANTEE



📅 29/08/2024 17:19:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)



## XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III



Ref. EKMWS07\_CIII

**Extra-long and low-slung scale for easy access to stretchers Medical Class III Version**

### Description

Thanks to its large dimensions and wide weighing range, is particularly suitable for weighing patients in wheelchairs or stretcher, as well as obese patients

- Large steel tray with non-slip and wear-resistant surface
- Low platform and integrated access ramps on both sides for easy access
- Stable and non-slip thanks to its height-adjustable rubber feet- Dynamic Weighing function for a stable weighing value, ideal for agitated patients

- Hold function: in case of service stoppage the weight display automatically "freezes" until the HOLD button is activated

- BMI function for calculating normal/overweight weight

- 2 casters and 2 handles make it easy to carry the scale and store it for a small

footprint- Analytical device with 2 magnets on the back, ideal for transport or fixing on metal surfaces

- The large LCD screen can be fixed to the wall thanks to the support supplied as a base or optional, placed on the Column EKMWS02 (H. +/- 95 cm)

- Reading: 0.100 kg

- Power supply: mains pack provided or batteries (6 x 1.5 V AA) not included, autonomy of about 50 h

- Dimensions: 1500 x, 860 x 68 mm (W x D x H)

- Weight: 40 kg

- Homologated in class III; Medical approval according to Directive 93/42/EEC

For any measuring device requiring medical Class III approval, the accessories must be ordered at the same time as the device so that the approval can be valid.

 29/08/2024 17:19:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)



## XXL WEIGHING PLATFORM (STRETCHER TROLLEY) KERN MWS-L 300KG-C MED III



Ref. **EKMWS07\_CIII**

### Features

CE Marking	CE0297	Medical device	Class IIa
Bar code	3661040173740	HS Code	84231090
Warranty (months)	36	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	86
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	150
Packed unit / Weight	55	Bare Unit / Weight (kg)	42
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 17:19:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)