

MWS-L WEIGHING PLATFORM

Ref. EKMWS07

KERN®



CE
0297

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 YEARS
GUARANTEE



29/08/2024 17:17:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MWS-L WEIGHING PLATFORM



Ref. EKMWS07

Extra-long and low-slung scale for easy stretcher access

Description

Thanks to its large dimensions and wide weighing range, is particularly suitable for weighing patients in wheelchairs or stretcher, as well as obese patients

- Large steel tray with non-slip and wear-resistant surface

- Low platform and integrated access ramps on both sides for easy access

- Stable and non-slip thanks to its height-adjustable

rubber feet- Dynamic Weighing function for a stable weighing value, ideal for agitated patients

-

Hold function: in case of service stoppage the weight display automatically "freezes" until the HOLD button is activated

- BMI function for calculating normal/overweight weight

- 2 casters and 2 handles make it easy to carry the scale and store it for a small

footprint- Analytical device with 2 magnets on the back, ideal for transport or fixing on metal surfaces

- The large LCD screen can be fixed to the wall thanks to the basic supplied bracket, or optionally, placed on the EKMWS02 Column (H. +/- 95 cm)-

Range: 300 kg- Reading:

0.100 kg

- Power supply: mains pack provided or batteries (6 x 1.5 V AA) not included, autonomy of about 50 h

-

Dimensions: 1500 x, 860 x 68 mm (W x D x H)

-

Weight: 40 kg

- Homologable in class III; Medical approval according to Directive 93/42/EEC NB: For a

version of this product approved Class III medical: reference EKMWS07_CIII

For any measuring device requiring a medical Class III approval, the accessories must be ordered at the same time as the device so that the approval can be valid.

29/08/2024 17:17:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MWS-L WEIGHING PLATFORM



Ref. EKMWS07

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040173733	HS Code	84231090
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	86
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	150
Packed unit / Weight	55	Bare Unit / Weight (kg)	42
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:17:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr