

BINAURALS FO MAGISTER II SILVER + TUBING BLACK (WITHOUT TIPS)



Ref. 506981



30/08/2024 18:16:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

**BINAURALS FO MAGISTER II SILVER +
TUBING BLACK (WITHOUT TIPS)**

Ref. 506981

Description

Tubing with lyre alu brush for stethoscope MAGISTER II from SPENGLER.
(without tips)

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030056	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.062
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:16:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr