

NON CHILL RING (BLACK) FOR BELL OF MAGISTER II SPENGLER



Ref. 507948

30/08/2024 18:52:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

NON CHILL RING (BLACK) FOR BELL OF MAGISTER II SPENGLER



Ref. 507948

Description

Spare part for stethoscopes MAGISTER® II from SPENGLER.
Black anti-cold ring for small bell.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030124	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:52:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr