

SINGLE TUBBING BLACK WITH BINAURAL ASSEMBLY FOR PEDIASONC SPENGLER



Ref. 507950



 30/08/2024 18:52:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SINGLE TUBBING BLACK WITH BINAURAL ASSEMBLY FOR PEDIASONC SPENGLER



Ref. 507950

Description

Spare part for Pediasonic® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446011277	HS Code	90181990
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.07
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:52:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr