

BLACK NYLON SOFT BAG SPENGLER, LARGE

Ref. 535019



04/11/2024 20:40:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BLACK NYLON SOFT BAG SPENGLER, LARGE



Ref. 535019

Description

Nylon pouch 25x20x5 cm with eclair® zip closure

Features

CE Marking	NONE	Medical device	No
Bar code	3700446002008	HS Code	90330000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.111
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

04/11/2024 20:40:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr