

SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS



Ref. 546011

CE
0123



30/08/2024 19:50:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS



Ref. 546011

Description

SPO2 pediatric sensor for SPENGLER spot-check station 2VS.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446029616	HS Code	90181990
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.056
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 19:50:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr