

EXTENSION TUBE FOR CUFFS FOR SPOT-CHECK 2VS




Ref. 546022



 30/08/2024 19:53:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

EXTENSION TUBE FOR CUFFS FOR SPOT-CHECK 2VS



Ref. 546022

Description

Extension cable cuffs for spot-check station 2VS from SPENGLER

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029517	HS Code	85049099
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.046
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:53:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr