

## CHARGEUR ALIMENTATION EU ET CORDON POUR SPOT-CHECK 2VS



Ref. 546023

30/08/2024 19:53:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CHARGEUR ALIMENTATION EU ET CORDON POUR SPOT-CHECK 2VS



Ref. 546023

### Description

European power charger and cord for SPENGLER spot-check station 2VS

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | NONE          | Medical device            | No       |
| Bar code                    | 3700446029692 | HS Code                   | 85049099 |
| Country                     | China         |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.093    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Packed unit / Quantity      | 1             |                           |          |

30/08/2024 19:53:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)