

## LIAN CUFF/1 TUBE/NYLON/BLACK/LARGE ADULT (L)



Ref. 333000



 30/08/2024 17:46:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN CUFF/1 TUBE/NYLON/BLACK/LARGE ADULT (L)



Ref. 333000

### Description

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color.

Washable at 30-40° C.

### Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | NONE          | Medical device            | Class I    |
| Bar code                    | 3700446008765 | HS Code                   | 9018901000 |
| Warranty (months)           | 6             | Country                   | Indonesia  |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.21       |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |

30/08/2024 17:46:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)