

NYLON CUFF/BLACK/1TUBE/CHILD S

Ref. 334820



30/08/2024 17:54:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

NYLON CUFF/BLACK/1TUBE/CHILD S



Ref. 334820

Description

Black Oxford nylon cuffs with quick coupling.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446006853	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.11
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:54:47 - For information purposes only, can be subject to change