

LIAN METAL / GREY COTTON CUFF / ADULT (M)

Ref. 521318



CE
0459



30/08/2024 19:37:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

Ambidextrous model offering great comfort of use thanks to its large dial.

LIAN METAL / GREY COTTON CUFF / ADULT (M)



Ref. 521318

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008840	HS Code	90189010
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.75
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:37:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr