

## PROBE FOR VASCULAR DOPPLER / 8 MHZ



Ref. 743312



30/08/2024 21:33:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PROBE FOR VASCULAR DOPPLER / 8 MHZ



Ref. 743312

### Description

8 MHz probe - cardiac probe for doppler

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446008161	HS Code	90181990
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.145
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:33:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)