

## DISPLAY OF 6 BLUE THERMOMETERS- TEMPO DUO

Ref. SP3100P



30/08/2024 21:37:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPLAY OF 6 BLUE THERMOMETERS-TEMPO DUO



Ref. SP3100P

The Tempo Duo allows a frontal and ear temperature measurement by infrared emissions

### Description

Display of 6 thermometers Tempo Duo

### Features

|                             |               |                           |           |
|-----------------------------|---------------|---------------------------|-----------|
| CE Marking                  | CE0123        | Medical device            | Class IIa |
| Bar code                    | 3700446023133 | HS Code                   | 90251920  |
| Warranty (months)           | 24            | Country                   | China     |
| DIMENSIONS AND WEIGHT       |               |                           |           |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc        |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 1.243     |
| SPECIFIC STORAGE CONDITIONS |               |                           |           |
| Minimum temperature         | NA            | Max temperature           | NA        |
| Minimum humidity            | NA            | Maximum humidity          | NA        |
| OTHER FEATURES              |               |                           |           |
| Packed unit / Quantity      | 1             |                           |           |

30/08/2024 21:37:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)