

## DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110P



30/08/2024 21:39:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110P

### Description

Cardboard counter display for 6 professional thermometers TEMPO PRO Spengler

Infrared technology allows the temperature to be taken on the forehead/temple without physical contact and effectively fights against cross-contamination.

With its "auto off" mode and "low battery indicator", no unpleasant surprises during your consultations.

TECHNICAL DATA• Measurement

mode: Non-contact infrared• Measuring distance: - 5 cm• Unit of measurement: C° / F°• Display range: 32.0°C – 43.0°C• Resolution: 0.1°C• Accuracy:

0.2°C•

Memory: 30 memories•

Measurement time: Less than 1 second• Screen:

LCD

OTHER

INFORMATION• Weight: 130 g•

Dimensions: 113(W) x 51 (W) x 162 (H)•

Automatic power saving power off

• Batteries supplied with (2 x AAA batteries)

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | NONE          | Medical device            | No       |
| Bar code                    | 3700446041939 | HS Code                   | 90251920 |
| Warranty (months)           | 24            | Country                   | China    |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | nc       |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Packed unit / Quantity      | 1             |                           |          |

30/08/2024 21:39:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)