

BATTERY FOR ES 60 / PRO M / PRO H SPENGLER



Ref. 550616



30/08/2024 19:57:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BATTERY FOR ES 60 / PRO M / PRO H SPENGLER



Ref. 550616

Description

SPENGLER ES 60 electronic tensiometer batteries
as well as SPENGLER's PRO M and PRO H

technology : Nickel-metal hydride

Voltage: 4.8V

Minimum capacity: 2500mAh

Features

CE Marking	NONE	Medical device	No
Bar code	3700446008901	HS Code	85075000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	5.1
Packed unit / Weight	0.12	Bare Unit / Weight (kg)	0.119
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:57:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr