

THERMAL PRINTER FOR 4VS SPENGLER MONITOR + PAPER ROLL



Ref. 551642



 30/08/2024 20:06:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

THERMAL PRINTER FOR 4VS SPENGLER MONITOR + PAPER ROLL



Ref. 551642

Description

SPENGLER 4VS Vital Signs Monitor Thermal Printer

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033118	HS Code	85049099
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.148
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:06:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr