

INFRARED THERMOMETER FOR 4VS MONITOR



Ref. 551644

 30/08/2024 20:07:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

INFRARED THERMOMETER FOR 4VS MONITOR



Ref. 551644

Description

Infrared Ear Replacement Thermometer for Vital Signs Monitor SPENGLER 4VS

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033132	HS Code	9025192090
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.162
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:07:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr